

A Technique for Characterizing Complex Polymer Solutions in Extensional Flows

Gavin Braithwaite
Stephen Spiegelberg

*Cambridge Polymer Group Inc.
11 Ward Street,
Somerville, MA 02143
<http://www.campoly.com>*



**Cambridge
Polymer Group** INC.

Consultation, Testing, and Instrumentation for Polymeric Materials



Introduction

- Industrial processing conditions always have some form of extensional deformation
 - Pumping, Filling, Spreading, Extrusion, Spinning and Ink-jets
- Standard characterization techniques all shear
- Some work on extensional deformation
 - Filament Stretching Extensional Rheometer (FiSER) and similar
 - Rheometrix RME
 - RFX Opposed Jet

V. Tirtaatmadja and T. Sridhar, *Journal of Rheology*, **1993**, 36, 3, 277-284

S. Spiegelberg, D. Ables, and G. McKinley, *Journal of Non-Newtonian Fluid Mechanics*, **1996**, 64, 2-3, 229-267



Motivation

- Clearly an experimental niche is open that requires an instrument that is:
 - compact (suitable for shop floor or lab bench operation)
 - robust (capable of fast turn around)
 - simple design (to allow easy operation and modification)
 - easy characterization (capable of rapidly yielding characteristic material parameters).
- Simplest intuitive test is “thumb -and-forefinger”



72nd SOR Hilton Head

3/16



Motivation

- Polymeric solutions common in industry
- In particular polymers derived from bio-materials common as process aids and rheological modifiers
- Materials such as those based on cellulose extremely common, finding uses in fields such as:
 - cosmetics (toothpastes, shampoos, lotions, adhesives,...)
 - foods (sauces, beverages,...)
 - pharmaceuticals (suspensions, creams, tablets,...)
 - paper/textiles (binders, thickeners,...)
 - ceramics (binders, thickeners,...)
- Hydroxy propyl methyl cellulose

72nd SOR Hilton Head

4/16

Capillary Breakup Extensional Rheometer (CaBER™)



- Original concept by Entov and co-workers
- Fluid sets time scale
 - “Instantaneous” stretch
 - Drainage governed by fluid
 - viscosity, elasticity, relaxation time, surface tension, effective extensional viscosity etc.
 - Gives measure of rheological parameters
 - Also measure of industrially relevant characteristic values
 - break-up time, stringiness etc.

References:


- A. V. Bazilevsky, V. M. Entov, and A. N. Rozhkov, “Liquid filament microrheometer and some of its applications,” presented at Proceedings of the 3rd European Rheology Conference, 1990.
- M. Renardy, *Journal of Non-Newtonian Fluid Mechanics*, **1994**, 51, 97-107
- D. Papageorgiou, *Physics of Fluids*, **1995**, 7, 7, 1529-1544
- V. M. Entov and E. J. Hinch, *J. Non-Newtonian Fluid Mech.*, **1997**, 72, 31-53
- A. V. Bazilevskii, V. M. Entov, M. M. Lerner, and A. N. Rozhkov, *Polym. Sci., Ser. A*, **1997**, 39, 3, 316-324
- J. Eggers, *Review of Modern Physics*, **1997**, 69, 3, 865-929
- V. M. Entov and E. J. Hinch, *J. Non-Newtonian Fluid Mech.*, **1997**, 72, 31-53
- A. Tripathi, P. Whittingstall, and G. H. McKinley, *Rheol. Acta*, **2000**, 39, 321-337
- G. H. McKinley and A. Tripathi, *Journal of Rheology*, 2000, 44, 3
- 72nd SOR Hilton Head

5/16

The CaBER™



- Laser Micrometer (resolution 5 μm , response time <1ms)
- Drive system:
 - Manual
 - Solenoid
 - Linear motor drive (stretch times from between 20 ms and >100s, linear or exponential)
- 6 mm diameter stainless steel plates
 - interchangeable
- Oven:
 - Low temperature to 120 °C
 - High temperature to 300 °C (in development)
- Force




72nd SOR Hilton Head

6/16

Test Fluids




- HPMC Methocel K100M from Dow Chemical
 - 0.5 wt% in water
 - Degree of substitution 1.4 (22% Methoxyl, 8.1% Hydroxypropyl)
 - approximate molecular weight 10^5 g/mol
 - Known to weakly gel under correct conditions
- PEO from Scientific Polymer Products
 - 0.5 wt% in water
 - molecular weight 2.5×10^6 g/mol
- Shear Rheology performed on TA Instruments AR1000N using a cone and plate geometry
- All experiments performed at $25^\circ\text{C} \pm 2^\circ\text{C}$



72nd SOR Hilton Head

7/16

Shear Rheology




72nd SOR Hilton Head

8/16



Shear Rheology




72nd SOR Hilton Head

9/16



Relaxation after steady shear




Data collected at 60 rad/s after a steady shear for 2 minutes of 2000 Pa

72nd SOR Hilton Head

10/16



Capillary Breakup




PEO in water

72nd SOR Hilton Head

11/16



Capillary Breakup




HPMC in water with and without Triton X114 surfactant (Union Carbide)

72nd SOR Hilton Head

12/16


Images



72nd SOR Hilton Head

13/16

Capillary breakup with time



Note: representative data show repeatability extremely difficult

72nd SOR Hilton Head

14/16

Summary



- PEO solution is a weakly elastic solution
 - CaBER relaxation time consistent with oscillatory data
 - breakup independent of stretch rate
- HPMC
 - oscillatory rheology resembles PEO data
 - shear rheology hints at more complexity (hysteresis, no plateau)
 - breakup appears to be rate sensitive
 - response time dependant
 - behavior strongly modified by surfactant
 - aggregates
 - structure
 - dehydration
 - Environment may be an issue

72nd SOR Hilton Head

15/16

<http://www.campoly.com>

Conclusions




- Demonstrated an instrument capable of studying capillary breakup of simple and complex fluids
- Showed comparison of polymeric solutions where polymers are interacting and associating
- Observed effects not obvious in shear rheology
- Clearly question mark over environmental conditioning of samples

72nd SOR Hilton Head

16/16



Evaporation



72nd SOR Hilton Head

17/16



CaBER™ - “Simple” Fluids

Constitutive Model	Form of Solution	Parameters found from regression to data
Newtonian, $\dot{\boldsymbol{\tau}} = \mathbf{h}, \mathbf{g}$	$D_{mid}(t) = 0.142 (s/h_0) (t_c - t)$	$t_c, s/h_0$
Power-Law Fluids $\dot{\boldsymbol{\tau}} = K \mathbf{g}^n$	$D_{mid}(t) = 2^{-n} (0.142) (s/K) (t_c - t)^n$	$t_c, s/K, n$
Upper Convected Maxwell $\dot{\boldsymbol{\tau}} + \lambda \dot{\boldsymbol{\tau}}_0 = \mathbf{h}, \mathbf{g}$	$D_{mid}(t) = D_0 (GD_0/s)^{1-n} \exp(-t/\lambda D_0)$	$t_c, G/s$

- Newtonian
 - dynamics of the drainage of fluid column and rupture of the liquid bridge governed by viscous and elastic properties of fluid

J. Eggers, *Review of Modern Physics*, **1997**, 69, 3, 865-929
 G. H. McKinley and A. Tripathi, *Journal of Rheology*, **2000**, 44, 3
 D. Papageorgiou, *Physics of Fluids*, **1995**, 7, 7, 1529-1544

- Viscoelasticity
 - rapid initial viscous-dominated phase, then intermediate time-scale where the filament drainage is governed by surface tension and elasticity

M. Renardy, *Journal of Non-Newtonian Fluid Mechanics*, **1994**, 51, 97-107
 V. M. Entov and E. J. Hinch, *J. Non-Newtonian Fluid Mech.*, **1997**, 72, 31-53

72nd SOR Hilton Head

18/16



CaBER™ - “Complex” Fluids

- Complex fluids
- A balance of forces on the fluid filament governs the evolution in the midpoint profile of the liquid bridge
- Allows derivation of an “apparent” extensional viscosity

$$3h_s \underbrace{\left(-\frac{2}{D_m} \frac{dD_m}{dt} \right)}_{\text{Viscous Stress}} = 3h \dot{\epsilon} = \underbrace{\frac{4F_z}{\pi D_m^2}}_{\text{Tensile Stress}} - \underbrace{\left[\frac{t_z - t_{rr}}{\rho} \right]}_{\text{Elastic/Non-Newtonian Stress}} - \underbrace{\frac{2s}{D_m}}_{\text{Capillary Pressure}}$$

$$\bar{\Pi}_{app}(\dot{\epsilon}) = \frac{2s/D_{mid}(t)}{\left\{ -\frac{2}{D_{mid}} \frac{dD_{mid}}{dt} \right\}} = \frac{s}{\frac{dD_{mid}}{dt}}$$


72nd SOR Hilton Head

19/16



Example Data


$$D_{mid}(t) = D_1 - \frac{(2X(t)-1)s}{3h_s} t$$



Glycerol

72nd SOR Hilton Head

$$D_{mid}(t) = \frac{GD_0}{s} \exp(-t/3I_c)$$




0.025 wt% $2.5 \times 10^6 \text{ Mw}$
 PS in Styrene oligomer

20/16



Example Images



Top Images:
Newtonian oil (styrene oligomer)

Bottom Images:
500 ppm Polystyrene (2.5×10^6 Mw)
in styrene oligomer

(From G. H. McKinley and
A. Tripathi, *J. Rheol.*, **2000**, 44, 3)

72nd SOR Hilton Head

21/16



Sample Images



72nd SOR Hilton Head

22/16