




Drug-Eluting Stents: Tracking Polymer Coating Degradation and API Release


Drug-eluting stents are designed to reduce restenosis, or re-narrowing of arteries after stent implantation. Typically, a polymer containing an active pharmaceutical ingredient (API) is coated on metallic stent struts, and the API is released over time to block cell growth within the implanted stent.

Beyond biocompatibility and safety assessment, drug-eluting stent characterization has historically focused on API release kinetics and stent/coating durability. Manufacturers typically demonstrate coating integrity by confirming the absence of cracks/peeling after expansion and fatigue. More recently, there has been increased emphasis on characterizing the fate of polymer excipients.

Cambridge Polymer Group (CPG) partnered with a cardiovascular company to understand the degradation behavior of a polymer coating on a drug-eluting stent.

In Vitro Degradation and Release

Coated stents were subjected to conditions designed to simulate physiological degradation and API release. The conditions considered variables such as temperature and conditioning media composition. Careful selection of aging conditions is critical to ensuring physiological relevance.



Changes in polymer coating with aging: (a) a broadening molar mass distribution and a reduction in the average molar mass and (b) a drop in the amount of polymer coating on the stents.

Gel Permeation Chromatography (GPC)


GPC is commonly utilized to characterize polymer molar mass distributions (MWD). Here, GPC was employed not only to evaluate changes in the molar mass distribution of the polymer coating but also to quantitate the amount of polymer coating and API remaining on the stent by exploiting refractive index and UV detection.

Results


GPC analysis showed a steady decrease in molar mass and broadening of the polymer MWD with increasing aging duration. Simultaneously, there was a loss of polymer and API mass on the stent, i.e., degradants were lost to the aging media and API was released from the coating. The API exhibited an initial burst release that transitioned to zero-order release with time. Preliminary results did not identify polymeric or oligomeric degradation products in the aging media, suggesting that the degradants may be small molecule compounds. Ongoing efforts are focused on characterizing the aging media using chromatography to provide a mass balance and identification of degradation products.

Key Points

- Understanding polymer degradation is critical to predicting long-term performance and safety of products over time.
- GPC is a powerful analytical tool for evaluating changes to polymer structure and concentration that occur under physiological conditions.



ANALYTICAL TESTING
BIOMEDICAL MATERIALS
MATERIALS CONSULTATION
RESEARCH & DEVELOPMENT



Cambridge Polymer Group, Inc. is a contract research laboratory specializing in materials. We partner with our clients to solve the world's toughest polymer problems utilizing our multi-disciplinary research team and full service laboratory.

We work with clients throughout the product life cycle to:

- Develop new materials
- Design prototypes for proof-of-concept studies
- Create and execute experimental design
- Validate and verify manufacturing processes
- Perform root-cause analysis in product failures

Cambridge Polymer Group, Inc. was founded in 1996 to provide a cost-effective resource for testing, research and development to clients who need periodic access to Ph.D.-level scientists and their support structure. We have developed a host of testing methods and materials for our clients, which number more than 1,000.

100 TradeCenter Drive, Suite 200, Woburn, Massachusetts 01801
P: 617-629-4400 • F: 617-629-9100 • info@campoly.com • www.campoly.com
ISO 17025 Accredited #3930.01 & ISO 9001 Certified #000912-1-US-1-QMS
DEA Licensed #RC0548606 & FDA Registered #3005793482

